

## COOL BAG FOR DIABETIQUE PM DIA, SILVER




Ref. MEDIA02\_B



 25/03/2025 23:22:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## COOL BAG FOR DIABETIQUE PM DIA, SILVER



Ref. MEDIA02\_B

Kit designed with the spectrum of patients with diabetes in mind

### Description

- Double insulated compartment to keep products cool
- Includes reusable cold gel
- Material: 420D fabric
- Weight: 130 g
- Exterior dimensions: 19 x 10 x 5 cm / 0.95 l
- Max. recommended capacity: 600 g
- Available in blue (ref. MEDIA02\_B) or silver (ref. MEDIA02\_A)
- Sold empty

### Features

CE Marking	NONE	Medical device	No
Bar code	8436564540690	HS Code	4202.12.99.90
Warranty (months)	24	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	12	Bare Unit / Width (cm)	10
Packed unit / Height (cm)	4	Bare Unit / Height (cm)	5
Packed unit / Length (cm)	20	Bare Unit / Length (cm)	19
Packed unit / Weight	0.214	Bare Unit / Weight (kg)	0.13
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Colors	Blue	Packed unit / Quantity	1

25/03/2025 23:22:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)