

LITESTICK PENLIGHT SPENGLER

Ref. 570620



25/03/2025 23:26:05 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

LITESTICK PENLIGHT SPENGLER



Ref. 570620

Description

This penlight is dedicated to general auscultation

- LED lighting for cold , intense and white light close to natural light
- Robust, made in aluminium
- light and compact
- on/off button
- 2 LR03 AAA batteries included

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446009830	HS Code	90189080907
Warranty (months)	24		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	5	Packed unit / Height (cm)	5
Packed unit / Length (cm)	20	Packed unit / Weight	0.046
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	10	Max temperature	40
OTHER FEATURES			
Colors	Black	Packed unit / Quantity	1

25/03/2025 23:26:05 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr