


## LITESTICK PENLIGHT SPENGLER


Ref. 570622



 25/03/2025 23:32:04 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## LITESTICK PENLIGHT SPENGLER



Ref. **570622**

### Description

This penlight is dedicated to general auscultation

- LED lighting for cold , intense and white light close to natural light
- Robust, made in aluminium
- light and compact
- on/off button
- 2 LR03 AAA batteries included

### Features

CE Marking	CE	Medical device	Class I
Bar code	3700446009878	HS Code	901890
Warranty (months)	24		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	5	Packed unit / Height (cm)	5
Packed unit / Length (cm)	20	Packed unit / Weight	0.046
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	10	Max temperature	40
<b>OTHER FEATURES</b>			
Colors	Green	Packed unit / Quantity	1

25/03/2025 23:32:04 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)