

6-CHANNEL ECG PAPER, ROLLED



Ref. 260025



18/03/2025 18:59:48 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

6-CHANNEL ECG PAPER, ROLLED



Ref. 260025

Description

Paper roll dimensions 110mm x 20m for ECG Cardiomate 6

Features

CE Marking	NONE	Medical device	No
Bar code	3700446027377	HS Code	4823400000
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	5	Packed unit / Height (cm)	5
Packed unit / Length (cm)	12	Packed unit / Weight	0.134
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

18/03/2025 18:59:48 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr