

MAGISTER STETHOSCOPE FOR CHILDREN / BLACK



Ref. 507210



05/03/2025 21:18:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

MAGISTER STETHOSCOPE FOR CHILDREN / BLACK



Ref. 507210

=====

Description

For pediatric auscultation. Stainless steel dual head chestpiece for pediatric use. Multifrequency diaphragm. Angled binaural with built-in spring and flexible, screw-on tips. Equipped with a non-chill ring for little ones' comfort.

Available in :

Blue : Ref 507214

Black: Ref 507210

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446009120	HS Code	90181990
Warranty (months)	60		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	16	Packed unit / Height (cm)	4
Packed unit / Length (cm)	32	Packed unit / Weight	0.375
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-10	Max temperature	40
Minimum humidity	NA	Maximum humidity	95
OTHER FEATURES			
Colors	Black	Packed unit / Quantity	1

📅 05/03/2025 21:18:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr