


## LUXASCOPE AURIS CCT LED 2.5 V - LUXAMED

LUXAMED®

Ref. A1.616.214



Technologie CCT  
2 couleurs de lumière



13/03/2025 23:23:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## LUXASCOPE AURIS CCT LED 2.5 V - LUXAMED

LUXAMED®

Ref. A1.616.214

Innovative otoscope in pocket size

### Description


- Innovative illumination: LED-Ring with 2x4 LEDs (concentric)  
switchable with one button press : warm white light or cold white light,  
no fibre optic,  
field of vision is not restricted
- Maintenance-free, because no lamp replacement necessary
  - Automatic switchoff after 3 minutes
  - Life expectancy of LEDs: approx. 100.000 hours due to ideal operation by constant current drive
  - Colour temperature: ring 1 approx. 3.000 K, ring 2 approx. 6.000 K
  - Optimized design of head
  - Pivoting 3-fold magnifying glass with locking function
  - Insufflation port for pneumatic test
  - Two-component handle (aluminium/fibre-glass reinforced), also combinable with Luxascope Dermatoscope CCT LED head
  - Electronic ON/OFF push-button switch
  - Power supply: 2x AAA Alkaline batteries, inclusive
  - Incl. 20 disposable ear funnels, grey, each 10 x 2.5 mm and 10 x 4.0 mm

Available in :  
Black : A1.616.114  
Blue : A1.616.214  
Grey : A1.616.314  
White : A1.616.914

 13/03/2025 23:23:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## LUXASCOPE AURIS CCT LED 2.5 V - LUXAMED

**LUXAMED®**

Ref. A1.616.214


### Features

CE Marking	CE	Medical device	Class I
Bar code	04260381963585	HS Code	90189084
Warranty (months)	24	Country	Germany
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	11	Packed unit / Height (cm)	4
Packed unit / Length (cm)	24	Packed unit / Weight	0.3
Bare Unit / Weight (kg)	0.12		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90
<b>OTHER FEATURES</b>			
Colors	Blue		

 13/03/2025 23:23:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)