

THE FRENCH TROLLEY FOR NC3/NC5/4VS/CT40


Ref. BDDSOCL01



 19/03/2025 21:42:38 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

THE FRENCH TROLLEY FOR NC3/NC5/4VS/CT40

Ref. **BDDSOCL01**

Very robust, extremely handy and "design"

Description

A French base made in France that helps us reduce our carbon footprint and demonstrates GSH's support in bringing manufacturing back to France!
Compatible with our monitors: NC3,

NC5, CT40, 4VS

Robust: designed in steel and weighing 13 kgs, with casters that can each withstand a load of 75kg.

Noble materials: steel + epoxy paint RAL 9016 for the foot and folded aluminum sheet for the basket

Extremely flexible mobility and maneuverability due to:
- A 5-pointed star equipped with 5 casters, 2 with brake.
- German-made castors (TENTE brand: EN 12530 standard)
- Nylon wheels with rubber tyre and precision metal ball bearing.

Features

CE Marking	NONE	Medical device	No
Bar code	3661040220130	HS Code	90189084
Country	France		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	30	Packed unit / Height (cm)	62
Bare Unit / Height (cm)	100	Packed unit / Length (cm)	100
Packed unit / Weight	17	Bare Unit / Weight (kg)	13
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 19/03/2025 21:42:38 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr