

## LUXASCOPE DERMATOSCOPE CCT LED 2.5 V - CONTACT PLATE WITH SCALE - LUXAMED ®

LUXAMED®


Ref. C1.716.114



13/03/2025 23:28:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## LUXASCOPE DERMATOSCOPE CCT LED 2.5 V - CONTACT PLATE WITH SCALE - LUXAMED®

**LUXAMED®**

Ref. **C1.716.114**

Innovative pocket-sized dermatoscope with CCT LED technology: choice of warm white or cool white light.

WITH graduated contact plate

### Description

- Lighting: LED ring with 2x5 LEDs (concentric), no fibre optics, field of vision is not limited.
- Long press on the switch to change from warm white to cool white
- Removable contact plate for easy cleaning
- Maintenance-free, as no bulb change is required.
- Continuous zoom with X10 magnification for a distortion-free image
- Optimised head design
- Automatic switch-off after 3 minutes
- LED life: approx. 100,000 hours at constant current
- Illuminance: approx. 10,000 Lux
- Colour temperature: ring 1: approx. 3000 K, ring 2: approx. 6000 K
- Handle made of two materials: aluminium and reinforced fibreglass, compatible with the head of the Luxascope Auris CCT LED otoscope
- ON/OFF push-button switch
- Power supply: 2x AAA alkaline batteries, included

Available in :

- Black : C1.716.114
- White: C1.716.914

### Features

CE Marking	NONE	Medical device	No
Bar code	04260381963622	HS Code	90189084
Warranty (months)	24	Country	Germany

For information purposes only, can be subject to change

#### DIMENSIONS AND WEIGHT

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

+33(0)4 42 90 31 31

info@gsh-med.fr

www.gsh-med.fr

OTHER FEATURES

