

Ref. E1.442.012



18/03/2025 19:20:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

Ref. E1.442.012

Description

McCOY blades for LUXAMED® laryngoscopes with integrated fibre optics are available in several sizes. The end of the optical fibre can be bent further.

	A (mm)	B(mm)
Size 2 (E1.442.012)	115	18
Size 3 (E1.443.012)	135	18
Size 4 (E1.444.012)	155	18

Laryngoscope blade with integrated fibre optics

Removable blade tip facilitates intubation and reduces the risk of injury.

Made from high-quality stainless steel for long life

Easy to clean and sterilise, as there are no dust traps

Maintenance-free

8,000 fibres for maximum life and light intensity

Optimum illumination with the largest possible cross-section (3 x 7mm)

Illuminance of around 14,000 Lux for 3.7V and 2.5V LEDs

Compatible with all handles conforming to ISO 7376 (eco-responsible standard)

Autoclavable up to 134°C (approx. 4,000 times)

Datamatrix with all important information for traceability

18/03/2025 19:20:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

Ref. E1.442.012

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381961871	HS Code	90189084
Warranty (months)	24		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90
OTHER FEATURES			
Size	Size 2		

18/03/2025 19:20:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr