


## GP DOCTOR BAG BLACK




Ref. MEGP02\_N



27/03/2025 21:28:28 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## GP DOCTOR BAG BLACK



Ref. **MEGP02\_N**

### Description

Lightweight and easy to use bag thanks to its internal velcro dividers that allow to organize the sanitary equipment in an easy and fast way. Its system of integrated shoulder straps makes it possible to transform it into a backpack.

- Light and resistant textile material, anti-tear (Ripstop)
- Special area for personal belongings, business cards, pen, etc ...
- Document compartment, A4
- Internal velcro dividers
- Insulated amp maker, 43 bulbs
- Zipper with eyelet for padlock.
- Large external pocket with double compartment and elastic bands
- Integrated shoulder straps for transformation into a backpack
- Weight: 1 kg
- Outdoor measures: 40 x 21 x 25 cm / 21 L
- Maximum recommended capacity: 3.5 Kg
- Sold empty
- Available in black (MEGP02\_N), anthracite/salmon (MEGP03\_R) and grey/turquoise (MEGP03\_B).

### Features

CE Marking	NONE	Medical device	No
Bar code	8436564540430	HS Code	4202129990
Warranty (months)	24	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	21
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	25
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	40
Packed unit / Weight	1.5	Bare Unit / Weight (kg)	1
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Colors	Black	Packed unit / Quantity	1

27/03/2025 21:28:28 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)