

Ref. M20500



DDM 

CE

DISPOSITIF MEDICAL
CLASSE
I
MEDICAL DEVICE

GARANTIE
2
2 YEARS
GUARANTEE

 28/03/2025 21:19:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr



=====

Ref. M20500

=====

Description

Pressure cuffs are medical devices designed for use in operating theatres, intensive care units and emergency transport vehicles. They are used to manually accelerate the infusion of blood and solutions when gravity is not sufficient to achieve the desired flow rates.

The GREY FUSE pressure cuff is compatible with all IV and blood bags, and comes in 3 sizes: 500cc (ref. M20500) , 1000cc (ref. M21000) and 3000cc (ref. M25000).

Optimised monitoring and reliability

Reusable and easy to clean/disinfect

Visualisation of the IV bag thanks to a transparent fabric made up of a screen and a polyamide mesh.

Pressure range from 0 to 300mmHg in 10mmHg increments

Pressure gauge accuracy ± 3 mmHg

Red zone on dial at 300mmHg to prevent overpressure

Pressure relief controlled by a simple tap.

All our cuff are made in our Castelculier factory near Agen.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700468211518	HS Code	981090
Warranty (months)	24	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	1.5
Bare Unit / Weight (kg)	380		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES 9 26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

+33(0)4 42 90 31 31

info@gsh-med.fr

www.gsh-med.fr