

SPHYGMOMANOMETER EASY 3 BLACK, BLADERLESS CUFF WITH VELCRO ADULT



Ref. DS99603_N

CE
0123

DISPOSITIF MEDICAL
CLASSE
Im
MEDICAL DEVICE

GARANTIE
1
YEAR
GUARANTEE



02/04/2025 01:00:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

SPHYGMOMANOMETER EASY 3 BLACK, BLADERLESS CUFF WITH VELCRO ADULT



Ref. DS99603_N

Anti-shock housing

Description

- Blood pressure monitor with cuff without pocket, without latex, washable nylon, black with velcro closure

- Graduation from 0 to 300 mmHg

- Shockproof housing

- Tube length: 65 cm

- Black rubber pear with ergonomic holding spoon

- Comes with nylon storage kit

- Cuff available in infant size (13 to 20 cm, ref DS99521_N) / child (18cm to 26 cm, ref DS99421_N) / adult (26 to 35 cm, ref DS99322_N) / strong arm (32cm to 42 cm, ref DS99603_N) / Multiarmband

Features

CE Marking	CE0123	Medical device	Class Im (measuring function)
Bar code	3661040173078	HS Code	90189020
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	18	Packed unit / Height (cm)	8
Packed unit / Length (cm)	12	Packed unit / Weight	0.5
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	10	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

02/04/2025 01:00:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr