

## DIAGNOSTIC PENLIGHT - WITH STANDARD BULB - LUXAMED

**LUXAMED**<sup>®</sup>

Ref. D1.212.512

CE



 04/04/2025 22:54:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## DIAGNOSTIC PENLIGHT - WITH STANDARD BULB - LUXAMED

**LUXAMED®**

**Ref. D1.212.512**

Penlight for eye and throat diagnostic., long life technology.

### Description

- Top quality ABS plastic housing
- Warm white light : comfortable for eyes
- Operating voltage 3.0 V
- Incl. 2 AAA batteries

Available in :

- White: D1.212.912
- Blue: D1.212.212
- Grey: D1.212.312
- Yellow: D1.212.512
- Kids : D1.312.912
- Black: D D1.212.112
- Red : D1.212.412
- Green : D1.212.612

### Features

CE Marking	CE	Medical device	Class I
Bar code	04260381961796	HS Code	90189084
Warranty (months)	24		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-10	Max temperature	40
Minimum humidity	10	Maximum humidity	90
<b>OTHER FEATURES</b>			
Colors	YELLOW		

📅 04/04/2025 22:54:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

🌐 [www.gsh-med.fr](http://www.gsh-med.fr)

☎ +33(0)4 42 90 31 31

✉ [info@gsh-med.fr](mailto:info@gsh-med.fr)